

2.5mm Fiber Optic Foam Swab

Product# 48040, 48042F

Product Description

Coventry™ Sealed Foam Swabs are constructed from reticulated 100 pores per inch medical grade foam which features the best open cell structure, providing excellent particle entrapment characteristics. The foam is thermally bonded to the swab handle without using adhesives. These swabs are the most economical and are ideal for general applications.

- Traps surface particles while cleaning
- High solvent capacity; holds solvent well
- Fiberless construction does not generate loose fibers or particles
- No adhesives or binders are used in the construction
- Economical for high quantity applications
- Sized to fit 2.5mm fiber optic connectors
- Soft polyurethane foam will not scratch surfaces

Typical Applications


Coventry™ Sealed Foam Swabs can be used to:

- Remove particulate and oily contaminants from 2.5mm fiber optic connector end faces
- Effective as end face cleaner in back plane applications when used in the Combination Cleaning Process with the Waterbased FiberWash Liquid and the QbE Wiper.
- To Clean excess adhesives after gluing
- Precision application of lubricants and solvents
- 2.5mm alignment sleeve cleaning
- Micro mechanical cleaning

Compatibility

Foam Heads are compatible with most common solvents such as isopropyl alcohol, and methanol. Not recommended for use with ketones such as acetone or methyl ethyl ketone.

Polypropylene Handles are generally compatible with all common solvents including most dilute or weak acids.




Availability

48040 2.75 " (6.8 cm) Polypropylene Flexible Tip
0.12" (3.0 mm) Foam Head, 500/ bag

48042F 2.75 " (6.8 cm) Polypropylene Flexible Tip
0.12" (3.0 mm) Foam Head, 50/ tube

2.5mm Fiber Optic Foam Swab

Product# 48040, 48042F



48040/48042F

Sealed Foam Swabs Test Data

Head Material	Sealed Foam
Swabs:	48040/48042F
Particle Counts:	
0.5 - 1 um	797
1.0 - 2 um	230
2.0 - 5 um	66
5.0 - 10 um	24
10 - 25 um	11
Over 25 um	2
TOTAL	1,130
NVR (ug/swab)	none detected
Chlorides (ug/swab)	none detected

The swabs were tested by Particle Measurement Technology. All the swabs tests were done using the recommended practices of the Institute of Environmental Science, Swabs and Wipes Working Group.

Technical and Application Assistance

Chemtronics provides a technical hotline to answer your technical and application related questions. The toll free number is: 1-800-TECH-401

Note:

This information is believed to be accurate. It is intended for professional end users having the skills to evaluate and use the data properly. CHEMTRONICS does not guarantee the accuracy of the data and assumes no liability in connection with damages incurred while using it.

Chemtronics® is a registered trademark of Chemtronics. All rights reserved. Coventry™ and Econowipes™ are trademarks of Chemtronics. All rights reserved.